


CD3D Antibody (R34179)

Catalog No.	Formulation	Size
R34179-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	915
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CD3D antibody is available for research use only.



Western blot testing of Jurkat [A] and MOLT4 [B] lysates with CD3D antibody at 1ug/ml.
Predicted molecular weight: ~19kDa.

Description

Additional name(s) for this target protein: CD3 delta

Application Notes

Optimal dilution of the CD3D antibody should be determined by the researcher.

Immunogen

Amino acids HETGRLSGAADTQ were used as the immunogen for this CD3D antibody. This sequence is specific to the alpha form of CD3 delta.

Storage

Aliquot and store the CD3D antibody at -20oC.