

## CD70 Antibody / CD70 molecule [clone CD70/12768R] (V6046)

Catalog No.	Formulation	Size
V6046-100UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	100 ug
V6046-20UG	0.2 mg/ml in 1X PBS with 0.05% BSA, 0.05% sodium azide	20 ug
V6046SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Recombinant Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG, kappa
<b>Clone Name</b>	CD70/12768R
<b>UniProt</b>	P32970
<b>Localization</b>	Cell membrane
<b>Applications</b>	Immunohistochemistry (FFPE) : 1-2ug/ml
<b>Limitations</b>	This CD70/CD70 molecule antibody is available for research use only.

### Application Notes

Optimal dilution of the CD70/CD70 molecule antibody should be determined by the researcher.

### Immunogen

A recombinant fragment corresponding to the C-terminus of human CD70 protein (exact sequence is proprietary) was used as the immunogen for the CD70/CD70 molecule antibody.

### Storage

CD70/CD70 molecule antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

