


CD83 Antibody (R36439)

Catalog No.	Formulation	Size
R36439-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9308
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This CD83 antibody is available for research use only.



Description

Additional name(s) for this target protein: B-cell activation protein; BL11; HB15

Application Notes

Optimal dilution of the CD83 antibody should be determined by the researcher.

Immunogen

Amino acids TTSCNSGTYR were used as the immunogen for this CD83 antibody. This sequence is common to isoforms

a/b/c.		
Storage Aliquot and store the CD83 antibody at -20	0oC.	