

CD8A Antibody / CD8 alpha chain [clone MSVA-008R] (V6143)

| Catalog No. | Formulation | Size |
|-------------|---|--------|
| V6143-100UG | Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide | 100 ug |
| V6143-20UG | Antibody in 1X PBS with 0.05% BSA, 0.05% sodium azide | 20 ug |

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

| | |
|---------------------------|--|
| Species Reactivity | Human |
| Format | Purified |
| Host | Rabbit |
| Clonality | Recombinant Rabbit Monoclonal |
| Isotype | Rabbit IgG, kappa |
| Clone Name | MSVA-008R |
| UniProt | P01732 |
| Localization | Cell membrane, Secreted |
| Applications | Immunohistochemistry (FFPE) : 1-2ug/ml |
| Limitations | This CD8A/CD8 alpha chain antibody is available for research use only. |

Application Notes

1. Optimal dilution of the CD8A/CD8 alpha chain antibody should be determined by the researcher.
2. This CD8A/CD8 alpha chain antibody is recombinantly produced by expression in human HEK293 cells.

Immunogen

A synthetic peptide corresponding to CD8alpha (within amino acids 135-235) was used as the immunogen for the CD8A/CD8 alpha chain antibody.

Storage

CD8A/CD8 alpha chain antibody with sodium azide - store at 2 to 8°C; antibody without sodium azide - store at -20 to -80°C.

