


## CD9 Antibody [clone CD9/7416] (V4260)


Catalog No.	Formulation	Size
V4260-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V4260-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V4260SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)


<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal (mouse origin)
<b>Isotype</b>	Mouse IgG1, kappa
<b>Clone Name</b>	CD9/7416
<b>Purity</b>	Protein A/G affinity
<b>UniProt</b>	P21926
<b>Localization</b>	Cell surface, Cytoplasm
<b>Applications</b>	Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT Western Blot : 2-4ug/ml
<b>Limitations</b>	This CD9 antibody is available for research use only.




IHC staining of FFPE human lymph node tissue with CD9 antibody (clone CD9/7416).  
HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



IHC staining of FFPE human pancreas tissue with CD9 antibody (clone CD9/7416). Inset: PBS used in place of primary Ab (secondary Ab negative control). HIER: boil tissue sections in pH 9 10mM Tris with 1mM EDTA for 20 min and allow to cool before testing.



SDS-PAGE analysis of purified, BSA-free CD9 antibody (clone CD9/7416) as confirmation of integrity and purity.



Western blot testing of human brain tissue lysate with CD9 antibody. Predicted molecular weight: 23-27 kDa.

## Description

CD9 is a type IV transmembrane glycoprotein with four transmembrane domains. CD9 on pre-B cells may play a role in cell-cell adhesion. In addition, CD9 may play a role in signal transduction mediated by interaction with low molecular weight GTP binding proteins. CD9 is expressed on early B cells, eosinophils, basophils and activated T cells and is a major component of the platelet cell surface. It is also expressed on most non-T acute lymphoblastic leukemia cells and on some acute myeloid and chronic lymphoid leukemia.

## Application Notes

Optimal dilution of the CD9 antibody should be determined by the researcher.

## Immunogen

Recombinant full-length human CD9 protein was used as the immunogen for the CD9 antibody.

## Storage

Aliquot the CD9 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.

