


CD90 Antibody / Thy1 (R34991)

Catalog No.	Formulation	Size
R34991-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Human, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7070
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CD90 antibody is available for research use only.



Western blot testing of human tonsil lysate with CD90 antibody at 0.1ug/ml. Predicted molecular weight 18~35kDa depending on glycosylation level.

Description

Additional name(s) for this target protein: THY1

Application Notes

Optimal dilution of the CD90 antibody should be determined by the researcher.

Immunogen

Amino acids PEHTYRSRTNFTSKY were used as the immunogen for this CD90 antibody.

Storage

Aliquot and store the CD90 antibody at -20oC.