


Cdc25A Antibody (F49547)

Catalog No.	Formulation	Size
F49547-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F49547-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Bovine, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P30304
Applications	Western Blot : 1:1000
Limitations	This Cdc25A antibody is available for research use only.



Cdc25A antibody western blot analysis in MCF-7 lysate. Predicted molecular weight: 59-70 kDa.

Description

CDC25A is a member of the CDC25 family of phosphatases. CDC25A is required for progression from G1 to the S phase of the cell cycle. It activates the cyclin-dependent kinase CDC2 by removing two phosphate groups. CDC25A is specifically degraded in response to DNA damage, which prevents cells with chromosomal abnormalities from progressing through cell division. CDC25A is an oncogene, although its exact role in oncogenesis has not been demonstrated.

Application Notes

Titration of the Cdc25A antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 57-84 from the human protein was used as the immunogen for this Cdc25A antibody.

Storage

Aliquot the Cdc25A antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.