


## CDC45 Antibody [clone AD-3] (RQ4490)

Catalog No.	Formulation	Size
RQ4490	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

[Bulk quote request](#)

<b>Availability</b>	1-2 weeks
<b>Species Reactivity</b>	Human
<b>Format</b>	Purified
<b>Host</b>	Rabbit
<b>Clonality</b>	Rabbit Monoclonal
<b>Isotype</b>	Rabbit IgG
<b>Clone Name</b>	AD-3
<b>Purity</b>	Protein A affinity
<b>UniProt</b>	O75419
<b>Localization</b>	Nucleus, Cytoplasm
<b>Applications</b>	Western Blot : 1:500-1:2000
<b>Limitations</b>	This CDC45 antibody is available for research use only.



Western blot testing of human Jurkat lysate with CDC45 antibody at 0.5ug/ml. Predicted molecular weight ~66 kDa.

### Description

Required for initiation of chromosomal DNA replication. [UniProt]

### Application Notes

Optimal dilution of the CDC45 antibody should be determined by the researcher.

## **Immunogen**

A synthetic peptide from human CDC45 was used as the immunogen for the CDC45 antibody antibody.

## **Storage**

Store the CDC45 antibody at -20oC.