


CDH11 Antibody (R33226)

Catalog No.	Formulation	Size
R33226-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Rat, Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1009
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This CDH11 antibody is available for research use only.



Western blot testing of HEK293 lysate overexpressing CDH11-MYC with CDH11 antibody [1ug/ml] in Lane A and with anti-MYC [1/1000] in lane C. Mock-transfected probed with anti-CDH11 [1ug/ml] in Lane B. Predicted molecular weight: ~88kDa but can be observed at 120-130kDa.

Description

Additional name(s) for this target protein: Cadherin 11; Osteoblast cadherin

Application Notes

Optimal dilution of the CDH11 antibody should be determined by the researcher.

Immunogen

Amino acids CESDQTKPLSNQP were used as the immunogen for this CDH11 antibody.

Storage

Aliquot and store the CDH11 antibody at -20oC.