


## CDK10 Antibody (R33124)

Catalog No.	Formulation	Size
R33124-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8558
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CDK10 antibody is available for research use only.



Western blot testing of human kidney lysate with CDK10 antibody at 2ug/ml. Predicted molecular weight: ~41kDa.

## Description

Additional name(s) for this target protein: Cyclin-dependent kinase 10

## Application Notes

Optimal dilution of the CDK10 antibody should be determined by the researcher.

## Immunogen

Amino acids AEPDLECEQIRLK were used as the immunogen for this CDK10 antibody.

## **Storage**

Aliquot and store the CDK10 antibody at -20oC.