


CDK10 Antibody (R33124)

Catalog No.	Formulation	Size
R33124-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8558
Applications	Western Blot : 2-4ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CDK10 antibody is available for research use only.



Western blot testing of human kidney lysate with CDK10 antibody at 2ug/ml. Predicted molecular weight: ~41kDa.

Description

Additional name(s) for this target protein: Cyclin-dependent kinase 10

Application Notes

Optimal dilution of the CDK10 antibody should be determined by the researcher.

Immunogen

Amino acids AEPDLECEQIRLK were used as the immunogen for this CDK10 antibody.

Storage

Aliquot and store the CDK10 antibody at -20oC.