


CDK7 Antibody (R34277)

Catalog No.	Formulation	Size
R34277-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1022
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CDK7 antibody is available for research use only.



Western blot testing of Jurkat lysate with CDK7 antibody at 0.3ug/ml. Predicted molecular weight: ~37kDa.

Description

Additional name(s) for this target protein: Cyclin-dependent kinase 7

Application Notes

Optimal dilution of the CDK7 antibody should be determined by the researcher.

Immunogen

Amino acids HRSEAKDGINRT were used as the immunogen for this CDK7 antibody.

Storage

Aliquot and store the CDK7 antibody at -20oC.