


## CDYL Antibody (R35977)

Catalog No.	Formulation	Size
R35977-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Pig, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	9425
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This CDYL antibody is available for research use only.



Western blot testing of human testis lysate with CDYL antibody at 1ug/ml. The expected ~75kDa band and the additional ~26kDa band are both blocked by the immunizing peptide.

## Description

Additional name(s) for this target protein: Chromodomain protein, Y-like

## Application Notes

Optimal dilution of the CDYL antibody should be determined by the researcher.

## Immunogen

Amino acids SMLKYLQRKIDEF were used as the immunogen for this CDYL antibody.

## **Storage**

Aliquot and store the CDYL antibody at -20oC.