


CEBPB Antibody (R33464)


Catalog No.	Formulation	Size
R33464-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat, Pig
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1051
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CEBPB antibody is available for research use only.



Western blot testing of mouse NIH3T3 lysate with CEBPB antibody at 1ug/ml. Predicted molecular weight: ~33/36-41 kDa (mouse/human).



Western blot testing of 1) human skin and 2) pig lung lysate with CEBPB antibody at 1ug/ml. Predicted molecular weight: 36-41 kDa.

Description

Additional name(s) for this target protein: CCAAT/enhancer binding protein (C/EBP) beta

Application Notes

Optimal dilution of the CEBPB antibody should be determined by the researcher.

Immunogen

Amino acids DHERAIDFSPYLEP were used as the immunogen for this CEBPB antibody.

Storage

Aliquot and store the CEBPB antibody at -20oC.