


Cebpz Antibody (R35776)

Catalog No.	Formulation	Size
R35776-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	12607
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Cebpz antibody is available for research use only.



Description

Additional name(s) for this target protein: CCAAT/enhancer binding protein zeta

Application Notes

Optimal dilution of the Cebpz antibody should be determined by the researcher.

Immunogen

Amino acids RSQLDDHPESDEEN were used as the immunogen for this Cebpz antibody. **Storage** Aliquot and store the Cebpz antibody at -20oC.