


## Cebpz Antibody (R35776)

Catalog No.	Formulation	Size
R35776-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	12607
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Cebpz antibody is available for research use only.



Western blot testing rat testis lysate with Cebpz antibody at 1ug/ml. Predicted molecular weight: ~120kDa.

## Description

Additional name(s) for this target protein: CCAAT/enhancer binding protein zeta

## Application Notes

Optimal dilution of the Cebpz antibody should be determined by the researcher.

## Immunogen

Amino acids RSQLDDHPESDEEN were used as the immunogen for this Cebpz antibody.

## Storage

Aliquot and store the Cebpz antibody at -20oC.