


CENTG1 Antibody (R35202)

Catalog No.	Formulation	Size
R35202-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q99490
Gene ID	116986
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This CENTG1 antibody is available for research use only.



Western blot testing of human amygdala lysate with CENTG1 antibody at 0.3ug/ml.
Predicted molecular weight: ~91kDa (PIKE-S/A isoform).

Description

Additional name(s) for this target protein: Centaurin gamma 1, ArfGAP with GTPase domain, ankyrin repeat and PH domain 2; AGAP2

Application Notes

Optimal dilution of the CENTG1 antibody should be determined by the researcher.

Immunogen

Amino acids QASLDSIREAVINSQ were used as the immunogen for this CENTG1 antibody. This sequence is specific to isoform 2/PIKE-S/PIKE-A.

Storage

Aliquot and store the CENTG1 antibody at -20oC.