


CFH Antibody (R34237)


Catalog No.	Formulation	Size
R34237-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	3075
Localization	Secreted
Applications	Western Blot : 0.03-0.1ug/ml IHC (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CFH antibody is available for research use only.



Western blot testing of human lung lysate with CFH antibody at 0.03ug/ml. Predicted molecular weight: ~139 kDa, observed here at ~170 kDa.



IHC testing of FFPE human liver tissue with CFH antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Complement factor H

Application Notes

Optimal dilution of the CFH antibody should be determined by the researcher.

1. This CFH antibody is specific for isoform a.

Immunogen

Amino acids HLPDRKKDQYK were used as the immunogen for this CFH antibody.

Storage

Aliquot and store the CFH antibody at -20oC.