


CHIT1 Antibody (R33796)

Catalog No.	Formulation	Size
R33796-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1118
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:2000
Limitations	This CHIT1 antibody is available for research use only.



Western blot testing of Daudi lysate with CHIT1 antibody at 1ug/ml. Predicted molecular weight: ~52kDa.

Description

Additional name(s) for this target protein: Chitinase 1

Application Notes

Optimal dilution of the CHIT1 antibody should be determined by the researcher.

Immunogen

Amino acids EHGPSPGQDTFCQ GK were used as the immunogen for this CHIT1 antibody.

Storage

Aliquot and store the CHIT1 antibody at -20oC.