


CHRM1 Antibody (R33511)

Catalog No.	Formulation	Size
R33511-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Mouse
Predicted Reactivity	Human, Rat, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1128
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CHRM1 antibody is available for research use only.



Description

Additional name(s) for this target protein: Cholinergic receptor, muscarinic 1; AChR

Application Notes

Optimal dilution of the CHRM1 antibody should be determined by the researcher.

Immunogen

Amino acids DKRRWRKIPKRPGSVH were used as the immunogen for this CHRM1 antibody. **Storage** Aliquot and store the CHRM1 antibody at -20oC.