


CHRNA5 Antibody (R34314)

Catalog No.	Formulation	Size
R34314-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1138
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CHRNA5 antibody is available for research use only.



Western blot testing of KELLY cell lysate with CHRNA5 antibody at 1ug/ml. Predicted molecular weight: ~53kDa.

Description

Additional name(s) for this target protein: Cholinergic receptor, nicotinic, alpha 5, LNCR2; CHRNA5

Application Notes

Optimal dilution of the CHRNA5 antibody should be determined by the researcher.

Immunogen

Amino acids HVDRYFTQKEET were used as the immunogen for this CHRNA5 antibody.

Storage

Aliquot and store the CHRNA5 antibody at -20oC.