


CHRNA5 Antibody (R34314)

Catalog No.	Formulation	Size
R34314-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1138
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CHRNA5 antibody is available for research use only.



Description

Additional name(s) for this target protein: Cholinergic receptor, nicotinic, alpha 5, LNCR2; CHRNA5

Application Notes

Optimal dilution of the CHRNA5 antibody should be determined by the researcher.

Immunogen

Amino acids HVDRYFTQKEET were used as the immunogen for this CHRNA5 antibody.

Storage Aliquot and store the CHRNA5 antibody at -20oC.	