


## CHRNA3 Antibody (R36044)

Catalog No.	Formulation	Size
R36044-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1142
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This CHRNA3 antibody is available for research use only.



Western blot testing of human cerebellum lysate with CHRNA3 antibody at 1ug/ml.  
Predicted molecular weight: ~53kDa.

## Description

Additional name(s) for this target protein: Cholinergic receptor, nicotinic, beta 3

## Application Notes

Optimal dilution of the CHRNA3 antibody should be determined by the researcher.

## Immunogen

Amino acids SRHVKKKEHFISQ were used as the immunogen for this CHRNA3 antibody.

## Storage

Aliquot and store the CHRNA3 antibody at -20oC.