


## CHRNA3 Antibody (R36044)

Catalog No.	Formulation	Size
R36044-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1142
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This CHRNA3 antibody is available for research use only.



Western blot testing of human cerebellum lysate with CHRNA3 antibody at 1 ug/ml.  
Predicted molecular weight: ~53kDa.

## Description

Additional name(s) for this target protein: Cholinergic receptor, nicotinic, beta 3

## Application Notes

Optimal dilution of the CHRNA3 antibody should be determined by the researcher.

## Immunogen

Amino acids SRHVKKEHFISQ were used as the immunogen for this CHRNA3 antibody.

## **Storage**

Aliquot and store the CHRNA3 antibody at -20oC.