


## CHUK Antibody (R33791)

Catalog No.	Formulation	Size
R33791-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1147
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CHUK antibody is available for research use only.



Western blot testing of HeLa nuclear lysate with CHUK antibody at 1ug/ml. Predicted molecular weight: ~85kDa.

## Description

Additional name(s) for this target protein: Conserved helix-loop-helix ubiquitous kinase

## Application Notes

Optimal dilution of the CHUK antibody should be determined by the researcher.

## Immunogen

Amino acids EHDHSLSCVVTPQD were used as the immunogen for this CHUK antibody.

## Storage

Aliquot and store the CHUK antibody at -20oC.