


CIRBP Antibody Biotin Conjugate (R35570BTN)

Catalog No.	Formulation	Size
R35570BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1153
Applications	Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:32000
Limitations	This CIRBP antibody is available for research use only.



Western blot testing of rat testis lysate with biotinylated CIRBP antibody at 1ug/ml.
Predicted molecular weight: ~18kDa.

Description

Additional name(s) for this target protein: Cold inducible RNA binding protein

Application Notes

Optimal dilution of the CIRBP antibody should be determined by the researcher.

Immunogen

Amino acids RDSYDSYATHNE were used as the immunogen for this CIRBP antibody.

Storage

Aliquot and store the CIRBP antibody at -20oC.