


CIRBP Antibody Biotin Conjugate (R35570BTN)

| Catalog No. | Formulation | Size |
|-------------|-----------------------------------------------------------------------------------------|--------|
| R35570BTN | $0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|----------------------|----------------------------------------------------------------------------------------------|
| Species Reactivity | Human, Mouse, Rat |
| Predicted Reactivity | Cow, Dog |
| Format | Biotin Conjugate |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 1153 |
| Applications | Western Blot : 1-3ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:32000 |
| Limitations | This CIRBP antibody is available for research use only. |



Description

Additional name(s) for this target protein: Cold inducible RNA binding protein

Application Notes

Optimal dilution of the CIRBP antibody should be determined by the researcher.

Immunogen

Amino acids RDSYDSYATHNE were used as the immunogen for this CIRBP antibody.

Storage

Aliquot and store the CIRBP antibody at -20oC.