


CIRBP Antibody (R35570)


Catalog No.	Formulation	Size
R35570-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1153
Localization	Cytoplasmic, nuclear
Applications	Western Blot : 0.01-0.03ug/ml Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CIRBP antibody is available for research use only.



Western blot testing of mouse testis lysate with CIRBP antibody at 0.01ug/ml. Predicted molecular weight: ~18 kDa.



IHC staining of FFPE human testis with CIRBP antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Cold inducible RNA binding protein

Application Notes

Optimal dilution of the CIRBP antibody should be determined by the researcher.

Immunogen

Amino acids RDSYDSYATHNE were used as the immunogen for this CIRBP antibody.

Storage

Aliquot and store the CIRBP antibody at -20oC.