


CIZ1 Antibody (R35393)

Catalog No.	Formulation	Size
R35393-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	25792
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CIZ1 antibody is available for research use only.



Western blot testing of HeLa nuclear lysate with CIZ1 antibody at 0.3ug/ml. Predicted molecular weight: ~91kDa.

Description

Additional name(s) for this target protein: Cip1-interacting zinc finger protein

Application Notes

Optimal dilution of the CIZ1 antibody should be determined by the researcher.

Immunogen

Amino acids QPNTQDKTPSKVTAR were used as the immunogen for this CIZ1 antibody.

Storage

Aliquot and store the CIZ1 antibody at -20oC.