


# CIZ1 Antibody (R35393)

Catalog No.	Formulation	Size
R35393-100UG	$0.5\ \text{mg/ml}$ in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	25792
Applications	Western Blot: 0.3-1ug/ml ELISA (peptide) LOD: 1:128000
Limitations	This CIZ1 antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Cip1-interacting zinc finger protein

#### **Application Notes**

Optimal dilution of the CIZ1 antibody should be determined by the researcher.

### **Immunogen**

Amino acids QPNTQDKTPSKVTAR were used as the immunogen for this CIZ1 antibody.

Storage Aliquot and store the CIZ1 antibody at -20oC.				