


Claudin 11 Antibody / CLDN11 [clone AEDI-3] (RQ4785)

Catalog No.	Formulation	Size
RQ4785	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Bulk quote request

Availability	1-2 weeks
Species Reactivity	Human, Mouse, Rat
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	AEDI-3
Purity	Affinity purified
UniProt	O75508
Applications	Western Blot : 1:500-1:2000
Limitations	This Claudin 11 antibody is available for research use only.



Description

Claudin 11 plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity. [UniProt]

Application Notes

Optimal dilution of the Claudin 11 antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human Claudin 11 / CLDN11 was used as the immunogen for the Claudin 11 antibody. **Storage** Store the Claudin 11 antibody at -20oC. Ordering: Phone: 858.663.9055 | Fax: 1.267.821.0800 | Email: info@nsjbio.com $\label{lem:copyright @NSJ Bioreagents.} \ \ \text{All rights reserved}.$