

Claudin 14 Antibody (R35844)

Catalog No.	Formulation	Size
R35844-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23562
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This Claudin 14 antibody is available for research use only.



Western blot testing of human liver lysate with Claudin 14 antibody at 0.3ug/ml.
Predicted molecular weight: ~25 kDa.

Description

Additional name(s) for this target protein: CLDN14

Application Notes

Optimal dilution of the Claudin 14 antibody should be determined by the researcher.

Immunogen

Amino acids SATHSGYRLNDYV were used as the immunogen for this Claudin 14 antibody.

Storage

Aliquot and store the Claudin 14 antibody at -20oC.