


CLEC1B Antibody (R34445)

Catalog No.	Formulation	Size
R34445-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	51266
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This CLEC1B antibody is available for research use only.



Western blot testing of human kidney lysate with CLEC1B antibody at 1ug/ml. Predicted molecular weight: ~27 kDa.

Description

Additional name(s) for this target protein: C-type lectin domain family 1, member B, CLEC2

Application Notes

Optimal dilution of the CLEC1B antibody should be determined by the researcher.

Immunogen

Amino acids KGHKCSPCDTNWRY were used as the immunogen for this CLEC1B antibody.

Storage

Aliquot and store the CLEC1B antibody at -20oC.