

CNBP Antibody (R34084)

Catalog No.	Formulation	Size
R34084-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7555
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This CNBP antibody is available for research use only.



Western blot testing of human muscle lysate with CNBP antibody at 2ug/ml. Predicted molecular weight: ~19kDa.

Description

Additional name(s) for this target protein: Cellular retroviral nucleic acid binding protein

Application Notes

Optimal dilution of the CNBP antibody should be determined by the researcher.

Immunogen

Amino acids GESGHLARECTIE were used as the immunogen for this CNBP antibody.

Storage

Aliquot and store the CNBP antibody at -20oC.