


## CNN2 Antibody (R33569)

Catalog No.	Formulation	Size
R33569-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1265
Applications	Western Blot : 0.03-0.1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This CNN2 antibody is available for research use only.



Western blot testing of HepG2 lysates with CNN2 antibody at 0.03ug/ml with [B] and without [A] blocking/immunizing peptide. Predicted molecular weight: ~34kDa.

## Description

Additional name(s) for this target protein: Calponin 2

## Application Notes

Optimal dilution of the CNN2 antibody should be determined by the researcher.

## Immunogen

Amino acids DPGEVPEYPPYYQEE were used as the immunogen for this CNN2 antibody.

## **Storage**

Aliquot and store the CNN2 antibody at -20oC.