


Cntn4 Antibody (R35744)

Catalog No.	Formulation	Size
R35744-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	152330
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Cntn4 antibody is available for research use only.



Western blot comparison of KO and WT tissue lysates with Cntn4 antibody at 0.5ug/ml. Lanes: KO mouse olfactory bulb (1) and cerebral cortex (2), and WT mouse olfactory bulb (3) and cerebral cortex (4). Lane 5: HEK293 lysate overexpressing mouse Cntn4.

Description

Additional name(s) for this target protein: Contactin 4

Application Notes

Optimal dilution of the Cntn4 antibody should be determined by the researcher.

Immunogen

Amino acids RRADGKPIARK were used as the immunogen for this Cntn4 antibody.

Storage

Aliquot and store the Cntn4 antibody at -20oC.