


## Coilin Antibody (R35680)

Catalog No.	Formulation	Size
R35680-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Rat
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8161
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Coilin antibody is available for research use only.



Western blot testing of MCF7 nuclear lysate with Coilin antibody at 1ug/ml. Predicted molecular weight: ~63kDa but routinely observed at 80kDa.

## Description

Additional name(s) for this target protein: CLN80, p80-coilin, COIL

## Application Notes

Optimal dilution of the Coilin antibody should be determined by the researcher.

## **Immunogen**

Amino acids RLVRDNDCLRVKLE were used as the immunogen for this Coilin antibody.

## **Storage**

Aliquot and store the Coilin antibody at -20oC.