


COL4A3 Antibody (R35621)

Catalog No.	Formulation	Size
R35621-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide $$	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1285
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This COL4A3 antibody is available for research use only.



IHC testing of FFPE human kidney with COL4A3 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Collagen, type IV, alpha 3

Application Notes

Optimal dilution of the COL4A3 antibody should be determined by the researcher.

Immunogen

Amino acids RGPTEYYDTYQEK were used as the immunogen for this COL4A3 antibody. **Storage** Aliquot and store the COL4A3 antibody at -20oC.