


COL4A3 Antibody (R35621)

Catalog No.	Formulation	Size
R35621-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1285
Applications	Immunohistochemistry (FFPE) : 4-6ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This COL4A3 antibody is available for research use only.



IHC testing of FFPE human kidney with COL4A3 antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Collagen, type IV, alpha 3

Application Notes

Optimal dilution of the COL4A3 antibody should be determined by the researcher.

Immunogen

Amino acids RGPTEYYDTYQEK were used as the immunogen for this COL4A3 antibody.

Storage

Aliquot and store the COL4A3 antibody at -20oC.