


Corneodesmosin Antibody (R35502)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35502-100UG | 0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

Bulk quote request

| Availability | 1-3 business days |
|--------------------|--|
| Species Reactivity | Human |
| Format | Antigen affinity purified |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 1041 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000 |
| Limitations | This Corneodesmosin antibody is available for research use only. |



Description

Additional name(s) for this target protein: CDSN

Application Notes

Optimal dilution of the Corneodesmosin antibody should be determined by the researcher.

Immunogen

Amino acids QVGNGSALPTNDNS were used as the immunogen for this Corneodesmosin antibody.

Storage Aliquot and store the Corneodesmosin antibody at -20oC.