


Corneodesmosin Antibody (R35502)

Catalog No.	Formulation	Size
R35502-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1041
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This Corneodesmosin antibody is available for research use only.



Western blot testing of cell line A431 lysate with Corneodesmosin antibody at 1ug/ml.
Predicted molecular weight: ~52kDa.

Description

Additional name(s) for this target protein: CDSN

Application Notes

Optimal dilution of the Corneodesmosin antibody should be determined by the researcher.

Immunogen

Amino acids QVGNGSALPTNDNS were used as the immunogen for this Corneodesmosin antibody.

Storage

Aliquot and store the Corneodesmosin antibody at -20oC.