


# **CORO1C Antibody (R35641)**

Catalog No.	Formulation	Size
R35641-100UG	0.5~mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

#### **Bulk quote request**

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23603
Applications	Western Blot : 0.2-1ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This CORO1C antibody is available for research use only.



## **Description**

Additional name(s) for this target protein: Coronin-3, Coronin 1C

## **Application Notes**

Optimal dilution of the CORO1C antibody should be determined by the researcher.

#### **Immunogen**

Amino acids RISKLEQQMAKIAA were used as the immunogen for this CORO1C antibody. **Storage** Aliquot and store the CORO1C antibody at -20oC.