


Cortactin Antibody (R36346)


Catalog No.	Formulation	Size
R36346-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat, Dog, Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	2017
Applications	Western Blot : 0.01-0.1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Cortactin antibody is available for research use only.



Western blot testing of mouse brain lysate with Cortactin antibody at 0.01ug/ml.
Molecular weight: routinely observed at ~80kDa.



Western blot testing of human 1) brain and 2) HEK293 lysate with Cortactin antibody at 0.1ug/ml and 0.03ug/ml, respectively. Molecular weight: routinely observed at ~80kDa.

Description

Additional name(s) for this target protein: CTTN; EMS1

Application Notes

Optimal dilution of the Cortactin antibody should be determined by the researcher.

Immunogen

Amino acids WKASAGHAVSIA were used as the immunogen for this Cortactin antibody.

Storage

Aliquot and store the Cortactin antibody at -20oC.