


## CPT1B Antibody (R33478)

Catalog No.	Formulation	Size
R33478-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

### Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1375
Localization	Cytoplasmic
Applications	Immunohistochemistry (FFPE) : 3-5ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This CPT1B antibody is available for research use only.



IHC testing of FFPE human heart with CPT1B antibody at 4ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

## Description

Additional name(s) for this target protein: Carnitine palmitoyltransferase 1B

## Application Notes

Optimal dilution of the CPT1B antibody should be determined by the researcher.

1. This CPT1B antibody will detect human isoforms a and c/d.

## **Immunogen**

Amino acids DIADLFQVPKAYS were used as the immunogen for this CPT1B antibody.

## **Storage**

Aliquot and store the CPT1B antibody at -20oC.