


CRB1 Antibody (RQ5239)


Catalog No.	Formulation	Size
RQ5239	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request


Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Host	Rabbit
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	P82279
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml Flow Cytometry : 1-3ug/million cells
Limitations	This CRB1 antibody is available for research use only.




Western blot testing of human 1) HEK293 and 2) U-87 MG cell lysate with CRB1 antibody. Predicted molecular weight ~154 kDa.




Flow cytometry testing of human U-87 MG cells with CRB1 antibody at 1ug/million cells (blocked with goat sera); Red=cells alone, Green=isotype control, Blue= CRB1 antibody.




IHC staining of FFPE human testis cancer with CRB1 antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 20 min and allow to cool before testing.



IHC staining of FFPE mouse testis with CRB1 antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 20 min and allow to cool before testing.



IHC staining of FFPE rat testis with CRB1 antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 20 min and allow to cool before testing.



IHC staining of FFPE human glioma with CRB1 antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 20 min and allow to cool before testing.

Description

Crumbs homolog 1 is a protein that in humans is encoded by the CRB1 gene. This gene encodes a protein which is similar to the *Drosophila* crumbs protein and localizes to the inner segment of mammalian photoreceptors. In *Drosophila* crumbs localizes to the stalk of the fly photoreceptor and may be a component of the molecular scaffold that controls proper development of polarity in the eye. Mutations in this gene are associated with a severe form of retinitis pigmentosa, RP12, and with Leber congenital amaurosis. Alternate splicing results in multiple transcript variants, some protein coding and some non-protein coding.

Application Notes

Optimal dilution of the CRB1 antibody should be determined by the researcher.

Immunogen

Amino acids FRTRDANVILHAEKEPEFLNISIQDSRLFFQLQ from the human protein were used as the immunogen for the CRB1 antibody.

Storage

Store the CRB1 antibody at -20oC.