


CRTC2 Antibody / TORC2 (R33565)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R33565-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|-----------------------------|--|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Mouse, Rat, Dog, Pig, Cow |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 200186 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:64000 |
| Limitations | This CRTC2 antibody is available for research use only. |



Western blot testing of Jurkat lysate with CRTC2 antibody at 1ug/ml. Predicted molecular weight: ~73 kDa, observed here at ~90 kDa.

Description

Additional name(s) for this target protein: CREB regulated transcription coactivator 2

Application Notes

Optimal dilution of the CRTC2 antibody should be determined by the researcher.

Immunogen

Amino acids DPAVEDSFRSDRLQ were used as the immunogen for this CRTC2 antibody.

Storage

Aliquot and store the CRTC2 antibody at -20oC.