


## CRTC2 Antibody / TORC2 (R33770)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33770-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                            |
|----------------------|------------------------------------------------------------|
| Availability         | 1-3 business days                                          |
| Species Reactivity   | Human                                                      |
| Predicted Reactivity | Mouse, Rat, Dog, Pig, Cow                                  |
| Format               | Antigen affinity purified                                  |
| Host                 | Goat                                                       |
| Clonality            | Polyclonal (goat origin)                                   |
| Isotype              | Goat Ig                                                    |
| Purity               | Antigen affinity                                           |
| Gene ID              | 200186                                                     |
| Applications         | Western Blot : 0.3-1ug/ml<br>ELISA (peptide) LOD : 1:32000 |
| Limitations          | This CRTC2 antibody is available for research use only.    |



Western blot testing of Jurkat lysate with CRTC2 antibody at 0.3ug/ml. Predicted molecular weight: ~73 kDa, observed here at ~85 kDa.

### Description

Additional name(s) for this target protein: TORC2, CREB regulated transcription coactivator 2

### Application Notes

Optimal dilution of the CRTC2 antibody should be determined by the researcher.

## Immunogen

Amino acids EFQSPLHSPLDSSR were used as the immunogen for this CRTC2 antibody.

## Storage

Aliquot and store the CRTC2 antibody at -20oC.