


CSF3R Antibody (R35094)

| Catalog No. | Formulation | Size |
|--------------|---|--------|
| R35094-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

| | |
|----------------------|---|
| Availability | 1-3 business days |
| Species Reactivity | Human |
| Predicted Reactivity | Cow |
| Format | Antigen affinity purified |
| Host | Goat |
| Clonality | Polyclonal (goat origin) |
| Isotype | Goat Ig |
| Purity | Antigen affinity |
| Gene ID | 1441 |
| Applications | Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000 |
| Limitations | This CSF3R antibody is available for research use only. |



Western blot testing of peripheral blood lymphocyte lysate with CSF3R antibody at 1ug/ml. Predicted molecular weight: ~92kDa.

Description

Additional name(s) for this target protein: Colony stimulating factor 3 receptor (granulocyte); CD114

Application Notes

Optimal dilution of the CSF3R antibody should be determined by the researcher.

Immunogen

Amino acids PPELGKSPLTHYT were used as the immunogen for this CSF3R antibody.

Storage

Aliquot and store the CSF3R antibody at -20oC.