


CSF3R Antibody (R35094)

Catalog No.	Formulation	Size
R35094-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow
Format	Antigen affinity purified
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1441
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This CSF3R antibody is available for research use only.



Western blot testing of peripheral blood lymphocyte lysate with CSF3R antibody at 1ug/ml. Predicted molecular weight: ~92kDa.

Description

Additional name(s) for this target protein: Colony stimulating factor 3 receptor (granulocyte); CD114

Application Notes

Optimal dilution of the CSF3R antibody should be determined by the researcher.

Immunogen

Amino acids PPELGKSPLTHYT were used as the immunogen for this CSF3R antibody.

Storage

Aliquot and store the CSF3R antibody at -20oC.