


CSK Antibody Biotin Conjugate (R33934BTN)

Catalog No.	Formulation	Size
R33934BTN	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog
Format	Biotin Conjugate
Host	Goat
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	1445
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : suitable ELISA (peptide) LOD : 1:128000
Limitations	This CSK antibody is available for research use only.



Western blot testing of Jurkat lysate with biotinylated CSK antibody at 0.5ug/ml.
Predicted molecular weight: ~50 kDa.

Description

Additional name(s) for this target protein: C-src tyrosine kinase

Application Notes

Optimal dilution of the CSK antibody should be determined by the researcher.

Immunogen

Amino acids EQLEHIKTHELHL were used as the immunogen for this CSK antibody.

Storage

Aliquot and store the CSK antibody at -20oC.